

# Oxygen Concentrators – Do's & Don'ts



## Do's:

- ❑ Use NEW mask/cannula and tubes for patient. Discard after Recovery. Ensure no bends in any tubing,
- ❑ Clean & sanitize oxygen concentrator & parts before use, follow clean up routine every 3 days.
- ❑ ONLY USE DISTILLED water for refilling the humidifier.
- ❑ Consult your doctor for queries on oxygen flow or if patient needs help upon using the oxygen concentrator.

## Dont's:

- ❑ Do not use Oxygen concentrator in closed spaces even if the AC is on. Place the oxygen concentrator away from walls and furniture in a well ventilated room.
- ❑ DO NOT USE concentrator without recommended break time.
- ❑ Do not aerosol products while using your portable oxygen concentrator.



**[www.oxygen4blr.in](http://www.oxygen4blr.in)**

For support/troubleshooting, Contact: **080 - 47187466** | **[info@oxygen4blr.in](mailto:info@oxygen4blr.in)**

Scan QR code for more information on how to use and maintain the concentrator

Web link of how to use an oxygen concentrator: <https://www.youtube.com/watch?v=mjcsf3JR92M>